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
N. Pearce, *et al.*, *Nat. Commun.* **2022**, *13*, 415.

Research Article
Rational Design and Assembly of Molecular Building Units into Metal-Organic Frameworks for Inverse Catalysis

Hakuhybotrol, a polyketide produced by Hypomyces pseudocorticicola, characterized with the assistance of 3D ED/MicroED
February 2023 · Organic & Biomolecular Chemistry 21(11)
DOI: 10.1039/d2ob02286a

Yoshihiro Watanabe · Shuhei Takahashi · Sho Ito · Masato Iwatsuki · Show all 15 authors · Synlett
DOI: 10.1055/s-0042-1751527

Electron Crystallography

Improving the Quality of Crystallographic Data

Liang Wang, M. Cristina D. Mallia, and D. ...

Communication

Fine Pore-Structure Engineering by Ligand Conformational Control of Naphthalene Diimide-Based Semiconducting Porous Coordination Polymers for Efficient Chemiresistive Gas Sensing

Dr. Ziqian Xue, Dr. Jia-Jia Zheng, Dr. Yusuke Nishiyama, Dr. Ming-Shui Yao, Dr. Yoshitaka Aoyama, Zeyu Fan, Dr. Ping Wang, Prof. Takashi Kajiwara, Prof. Yoshiki Kubo, and Prof. ...
First published: 15 November 2022 | <https://doi.org/10.1002/anie.202211111>

Aljaž Škrjanc, Martin Gazvoda, and Nataša Zupan

White qbe Rod Se...
Chemical Synthesis

research papers

A Journal of the German Chemical Society

International Porous Ladder...
Wang, and Zhenjie Zhang*

Catalysis in Germany

Alternative Epoxidation Reactions Catalyzed by...
Complexes with a Hydrogen-Bonding Site

An oxidative photocyclization approach to the synthesis of Securiflustra securifrons alkaloids

BRANDON W. ALEXANDER · NOAH M. BARTFIELD · VAANI GUPTA · BRANDON Q. MERCADO · MARK DEL CAMPO · AND SETH B. HERZON

SCIENCE · 22 Feb 2024 · Vol 383, Issue 6685 · pp. 849-854 · DOI: 10.1126/science.adg167

Christian R. Göb...

ROYAL SOCIETY OF CHEMISTRY

View Article Online
View Journal | View Issue

Spin-state dependent pressure responsiveness of Fe(II)-based triazolate metal-organic frameworks*

Silva M. Kronawitter, Richard Röß-Ohlenroth, Sebastian A. Hallweger, Marcel Hirle, Hans-Albrecht Krug von Nidda, Tobias Luxenhofer, Emily Myatt, Jem Pitcairn, Matthew J. Cliffe, Dominik Daisenberger, Jakob Wojciechowski, Dirk Volkmer, and Gregor Kieslich

Enna...
Abstraction and...

Yoshitaka Ashitani, Yoshiki Kubo, Ken-ichi Otake, and Ken-ichi Otake

GDCh

From a Molecular Complex to a M-HOF and a MOF

Andrichs, and Norbert Stock*

6/6
8
Citations
LEARN ABOUT THESE METRICS

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