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The screenshot shows the 'Cluster Analysis Browser' window with a dendrogram on the left. The 'Base Profile' tab is selected. Below the dendrogram, there are several settings panels:

- Conditions:** Conversion: Data reduction, Data reduction condition: Edit..., Correlation type: Default, Range: Min: 2,0000, Max: 60,0000, Calculate.
- View settings:** Show tooltip (checked), Original (selected), Converted, Show full path.
- Base decomposition settings:** Component: 3, Iteration: 1000, Calculate, EditDataset: Edit...

Red boxes highlight the 'Conditions' and 'Base decomposition settings' panels. A red arrow points from the 'Calculate' button in the 'Base decomposition settings' panel to the 'EditDataset' button.

The screenshot shows the 'Cluster Analysis Browser' window with a 'Base Profile' plot. The plot shows Intensity (cps) vs. 2θ (ω, °) with multiple peaks. The 'Base Profile' tab is selected. Below the plot, there are several settings panels:

- Conditions:** Conversion: Original Data, Data reduction condition: Edit..., Correlation type: Default, Range: Min: 2,0000, Max: 60,0000, Calculate.
- View settings:** Show tooltip (checked), Original (selected), Converted, Show full path.
- Base decomposition settings:** Method: Semi-supervised minimization, Component: 3, Estimate, Iteration: 1000, Calculate, EditDataset: Edit...

Red boxes highlight the 'Base decomposition settings' panel. A red arrow points from the 'Calculate' button in the 'Base decomposition settings' panel to the 'EditDataset' button.

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