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Attachment: ASC Robot

Measurement conditions

Suffix measurement date and time to filename  Suffix measurement part type name to filename  Always apply sample name to file name  Add file name defined in each part activity to file name

Copy to All Cassette Copy to Same Cassette Copy to Selected Sample

Import... Export... Options... Clear All Cassette Clear Selected Cassette

No.	Status	Exec.	Interrupt	Flow for Measurement	Sample Name	Folder	File Name	Memo
▶ Cassette #1 Interrupt order: 1 There are no samples to be measured.								
▲ Cassette #2 Interrupt order: 2 There are 5 samples to be measured.								
1	Standby	<input checked="" type="checkbox"/>	<input type="checkbox"/>	C:\Users\User1\Documents\%a-sasaki\...	ABP21B001-1	C:\Users\User1\Documents\%a-sasaki\...	ABP21B001-1-2-01.rasx	
2	Standby	<input checked="" type="checkbox"/>	<input type="checkbox"/>	C:\Users\User1\Documents\%a-sasaki\...	ABP21B001-2	C:\Users\User1\Documents\%a-sasaki\...	ABP21B001-2-2-02.rasx	
3	Standby	<input checked="" type="checkbox"/>	<input type="checkbox"/>	C:\Users\User1\Documents\%a-sasaki\...	ABP21B001-3	C:\Users\User1\Documents\%a-sasaki\...	ABP21B001-3-2-03.rasx	
4	Standby	<input checked="" type="checkbox"/>	<input type="checkbox"/>	C:\Users\User1\Documents\%a-sasaki\...	ABP21B001-4	C:\Users\User1\Documents\%a-sasaki\...	ABP21B001-4-2-04.rasx	
5	Standby	<input checked="" type="checkbox"/>	<input type="checkbox"/>	C:\Users\User1\Documents\%a-sasaki\...	ABP21B001-5	C:\Users\User1\Documents\%a-sasaki\...	ABP21B001-5-2-05.rasx	
6	Skipped	<input type="checkbox"/>	<input type="checkbox"/>					
7	Skipped	<input type="checkbox"/>	<input type="checkbox"/>					
8	Skipped	<input type="checkbox"/>	<input type="checkbox"/>					
9	Skipped	<input type="checkbox"/>	<input type="checkbox"/>					
10	Skipped	<input type="checkbox"/>	<input type="checkbox"/>					
11	Skipped	<input type="checkbox"/>	<input type="checkbox"/>					
12	Skipped	<input type="checkbox"/>	<input type="checkbox"/>					
13	Skipped	<input type="checkbox"/>	<input type="checkbox"/>					
14	Skipped	<input type="checkbox"/>	<input type="checkbox"/>					
15	Skipped	<input type="checkbox"/>	<input type="checkbox"/>					
16	Skipped	<input type="checkbox"/>	<input type="checkbox"/>					
17	Skipped	<input type="checkbox"/>	<input type="checkbox"/>					

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