

Appendix A. Description of Goniometers and Aspects of Data Collection

Tables A.I and A.II list the details of the experimental conditions for the X-ray measurements. Briefly, a 4-axis (Φ , χ , Ω , 2θ) goniometer was employed for the stress determinations using both the Ω - and Ψ -goniometer geometries [15]. During scanning, the specimens were oscillated ± 2 mm in plane to improve particle statistics. No stress determinations were performed on the 2-axis (θ - θ) goniometer.

Specimen alignment was accomplished using a dial gauge probe, which was accurate to $\pm 5 \mu\text{m}$ and a telescope. Here, the relative dis-

tance to the center of rotation is known, and the diffracting surface is positioned accordingly. Further, a telescope was initially employed in specimen alignment. The position of the specimen was confirmed by rotating 180° about an axis parallel to the diffracting surface of the specimen and observing that this surface was coincident with the horizontal cross hair at both -90 and $+90^\circ \chi$. Goniometer alignment was ensured by examining LaB_6 powder on a zero background plate. The maximum observed peak shift for the (510) reflection of LaB_6 ($141.7^\circ 2\theta$) was less than $0.03^\circ 2\theta$ for Ω and χ tilting as described in Table A.I.

Table A.I. Experimental conditions of the X-ray measurements 4-axis unit.*

| Parameter | Condition |
|--|---|
| Equipment | PTS goniometer Spellman DF3 series 4.0 kW generator Liquid N ₂ -cooled Ge detector |
| Power | 1.32 kW; 40 kV, 33 mA (long fine focus 0.4x12mm) |
| Radiation | Cu, $\lambda = 1.54056 \text{ \AA}$ |
| Incidence slit divergence | 0.12° |
| Receiving slit acceptance | 0.25° ; radial divergence limiting (RDL) Soller slit |
| Source to specimen distance | 290 mm |
| Specimen to back slit distance | 290 mm |
| Tilt axis and angles | χ & Ω ; $0, \pm 55^\circ$ |
| Scans | 0.02 & $0.05^\circ 2\theta/\text{step}$ |
| *Software: DMSNT v1.39-1B (build 125). | |

Table A.II. Experimental conditions of the X-ray measurements 2-axis unit.**

| Parameter | Condition |
|---|---|
| Equipment | θ - θ goniometer Seifert ID-3000, 3.5kW generator Bicron scintillation detector with curved graphite diffracted beam monochromator*** |
| Power | 1.8 kW; 45 kV, 40 mA (normal focus 1 x 10 mm) |
| Radiation | Cu, $\lambda = 1.54056 \text{ \AA}$ |
| Incidence slit divergence | 0.7° |
| Receiving slit acceptance | 0.3 mm |
| Source to specimen distance | 250 mm |
| Specimen to back slit distance | 250 mm |
| Scans | $0.01^\circ 2\theta/\text{step}$; $1-2^\circ/\text{min}$. |
| **Software: DMSNT v1.37. | |
| ***Advanced Ceramics Corp., pyrolytic graphite, ZYA grade, (002) oriented, $0.4^\circ \pm 0.1^\circ$ mosaic spread, $r = 225$ mm. | |



Fig. A1. 4-axis (base θ - 2θ) and 2-axis (θ - θ) goniometers.



Fig. A2. 4-axis unit: Goniometer movements for χ tilting and symmetric diffraction ($\Omega=2\theta/2$ always). Left to right χ : $-55, 0, +55^\circ$.



Fig. A3. 4-axis unit: Goniometer movements for Ω tilting and asymmetric diffraction ($\Omega \neq 2\theta/2$). Left to right Ω : $25.9, 70.9, 115.9^\circ$ corresponding to Ψ : $-55, 0, +55^\circ$.




Fig. A4. Sample mounts for the (A) 4-axis and (B) $\theta-\theta$ units. The 4-axis unit uses a dial gage probe with $\pm 5 \mu\text{m}$ precision while the $\theta-\theta$ unit uses a fiducial surface with no easy adjustments.




Fig. A5. $\theta-\theta$ unit: Photo showing the utilization of a laser pointer to find the center of rotation of the goniometer.

Appendix B. Diffraction Angle Calibration

PURPOSE

The Diffraction Angle Calibration provides a means to correct the angular position of the diffraction lines for both instrumental and physical aberrations. This correction is useful when accurate peak positions are to be determined, e.g., calculating lattice parameters, as well as improving phase identification results.

CALIBRATION STEPS

1. Prepare a slurry mount of a standard/reference powder ~ 325 mesh particle size (e.g., LaB_6 , Si, Al_2O_3) with either methanol or acetone and spread on a sample support, preferably a zero background plate. An Al_2O_3 or quartz plate each of random/near random texture may also be used. For the high temperature unit, the support will be a properly mounted heater strip.
2. Select the divergence and receiving slits typical of those used for samples to be measured on the unit. The irradiated sample length

Table B.I. Lattice parameter refinement of LaB₆ data in Fig. B1 taken on the 4-axis unit.

| [D07893.raw] Lab6 on zero bkgd plate | | | | | | | | | Cell Refinement Report | |
|--|---------|---------|---------|---------|--------|---------|---------|---------|------------------------|-------|
| SCAN: 15.0/155.0/0.02/0.2(sec), Cu, I(max)=5129, 06/21/01 08:30 | | | | | | | | | | |
| PEAK: 25-pts/Parabolic Filter, Threshold=1.0, Cutoff=0.05%, BG=3/1.0, Peak-Top=Summit | | | | | | | | | | |
| NOTE: Intensity = Counts, 2T(0)=0.0(°), Wavelength to Compute d-Spacing = 1.54056A (Cu/K-alpha1) | | | | | | | | | | |
| Cell Type = Cubic, Pm3m (221) | | | | | | | | | | |
| Refined Cell = 4.15984(0.000402) [] | | | | | | | | | | |
| Vol= 71.98 A^3, Density(c)= 4.7002 (Chemical Formula =LaB6, Z=1.0) | | | | | | | | | | |
| 2-Theta Error Window = 0.3(°), Zero Offset = 0.0(°), Displacement = 0.0(°) | | | | | | | | | | |
| ESD of Fit = 0.0882(°), lDelta 2-Theta = 0.0749(°), lDelta d = 0.00317(A), F(24) = 13.4(24) | | | | | | | | | | |
| Intensity Weighting = Sqrt(I%), 2-Theta Range = 15.0/155.0(°), (x) Outlier Rejection at 2.0 sigmas | | | | | | | | | | |
| 2T(cor) = 2T(obs) - Zero Offset - Displacement (cal=Calculated, obs=Observed, cor=Corrected) | | | | | | | | | | |
| # | (h k l) | 2T(cal) | 2T(cor) | 2T(obs) | Delta | d(cal) | d(cor) | d(obs) | Del-d | I% |
| 1 | (1 0 0) | 21.342 | 21.185 | 21.185 | 0.157 | 4.15984 | 4.19028 | 4.19028 | -0.03043 | 56.9 |
| 2 | (1 1 0) | 30.362 | 30.209 | 30.209 | 0.153 | 2.94145 | 2.95602 | 2.95602 | -0.01457 | 97.7 |
| 3 | (1 1 1) | 37.414 | 37.278 | 37.278 | 0.136 | 2.40168 | 2.41013 | 2.41013 | -0.00844 | 47.0 |
| 4 | (2 0 0) | 43.473 | 43.354 | 43.354 | 0.119 | 2.07992 | 2.08537 | 2.08537 | -0.00545 | 25.4 |
| 5 | (2 1 0) | 48.920 | 48.809 | 48.809 | 0.111 | 1.86034 | 1.86430 | 1.86430 | -0.00397 | 49.0 |
| 6 | (2 1 1) | 53.946 | 53.844 | 53.844 | 0.102 | 1.69825 | 1.70122 | 1.70122 | -0.00298 | 30.8 |
| 7 | (2 2 0) | 63.167 | 63.059 | 63.059 | 0.108 | 1.47073 | 1.47298 | 1.47298 | -0.00226 | 10.3 |
| 8 | (2 2 1) | 67.492 | 67.388 | 67.388 | 0.104 | 1.38661 | 1.38849 | 1.38849 | -0.00188 | 28.6 |
| 9 | (3 1 0) | 71.685 | 71.593 | 71.593 | 0.093 | 1.31546 | 1.31693 | 1.31693 | -0.00147 | 19.5 |
| 10 | (3 1 1) | 75.779 | 75.702 | 75.702 | 0.077 | 1.25424 | 1.25532 | 1.25532 | -0.00108 | 14.3 |
| 11 | (2 2 2) | 79.800 | 79.720 | 79.720 | 0.080 | 1.20084 | 1.20184 | 1.20184 | -0.00100 | 2.3 |
| 12 | (3 2 0) | 83.771 | 83.694 | 83.694 | 0.077 | 1.15373 | 1.15460 | 1.15460 | -0.00087 | 9.9 |
| 13 | (3 2 1) | 87.711 | 87.717 | 87.717 | -0.005 | 1.11176 | 1.11171 | 1.11171 | 0.00005 | 100.0 |
| 14 | (4 0 0) | 95.579 | 95.514 | 95.514 | 0.065 | 1.03996 | 1.04050 | 1.04050 | -0.00054 | 4.2 |
| 15 | (3 2 2) | 99.544 | 99.490 | 99.490 | 0.054 | 1.00891 | 1.00931 | 1.00931 | -0.00040 | 19.1 |
| 16 | (3 3 0) | 103.555 | 103.509 | 103.509 | 0.045 | 0.98048 | 0.98079 | 0.98079 | -0.00031 | 15.0 |
| 17 | (3 3 1) | 107.635 | 107.603 | 107.603 | 0.032 | 0.95433 | 0.95453 | 0.95453 | -0.00019 | 7.0 |
| 18 | (4 2 0) | 111.810 | 111.784 | 111.784 | 0.026 | 0.93017 | 0.93031 | 0.93031 | -0.00014 | 10.5 |
| 19 | (4 2 1) | 116.110 | 116.087 | 116.087 | 0.023 | 0.90775 | 0.90787 | 0.90787 | -0.00011 | 18.4 |
| 20 | (3 3 2) | 120.576 | 120.578 | 120.578 | -0.002 | 0.88688 | 0.88687 | 0.88687 | 0.00001 | 10.2 |
| 21 | (4 2 2) | 130.228 | 130.267 | 130.267 | -0.039 | 0.84912 | 0.84899 | 0.84899 | 0.00014 | 8.4 |
| 22 | (4 3 0) | 135.595 | 135.654 | 135.654 | -0.059 | 0.83197 | 0.83179 | 0.83179 | 0.00018 | 8.9 |
| 23 | (5 1 0) | 141.535 | 141.625 | 141.625 | -0.090 | 0.81581 | 0.81559 | 0.81559 | 0.00022 | 36.2 |
| 24 | (3 3 3) | 148.382 | 148.421 | 148.421 | -0.039 | 0.80056 | 0.80048 | 0.80048 | 0.00008 | 11.9 |




Fig. B1. 4-axis unit: Diffraction pattern of LaB₆ shows peak positions are systematically low relative to PDF card #34-427.

Table B.II. Lattice parameter refinement of LaB₆ data in Fig. B2 taken on the θ - θ unit.

| [X01593.raw] LaB6 std., Sollers in place | | | | | | | Cell Refinement Report | | | |
|---|---------|---------|---------|---------|--------|---------|------------------------|---------|------------|-------|
| SCAN: 10.0/159.99/0.01/0.2(sec), Cu, I(max)=1866, 06/22/01 15:24 | | | | | | | | | | |
| PEAK: 21-pts/Parabolic Filter, Threshold=1.0, Cutoff=0.05%, BG=3/1.0, Peak-Top=Summit | | | | | | | | | | |
| NOTE: Intensity = Counts, 2T(0)=0.0(°), Wavelength to Compute d-Spacing = 1.54056A (Cu/K-alpha1) | | | | | | | | | | |
| Cell Type = Cubic, Pm3m (221) | | | | | | | | | | |
| Refined Cell = 4.15179(0.000177) Å | | | | | | | | | | |
| Vol= 71.57 Å ³ , Density(c)= 4.7276 (Chemical Formula =LaB6, Z=1.0) | | | | | | | | | | |
| 2-Theta Error Window = 0.3(°), Zero Offset = 0.0(°), Displacement = 0.0(°) | | | | | | | | | | |
| ESD of Fit = 0.0118(°), IDelta 2-Theta = 0.0058(°), IDelta dI = 0.00014(Å), F(16) = 171.0(16) | | | | | | | | | | |
| Intensity Weighting = Sqrt(I%), 2-Theta Range = 10.0/159.99(°), (x) Outlier Rejection at 2.0 sigmas | | | | | | | | | | |
| 2T(cor) = 2T(obs) - Zero Offset - Displacement (cal=Calculated, obs=Observed, cor=Corrected) | | | | | | | | | | |
| # | (h k l) | 2T(cal) | 2T(cor) | 2T(obs) | Delta | d(cal) | d(cor) | d(obs) | Del-d | % |
| 1 | (1 0 0) | 21.384 | 21.384 | 21.384 | 0.000 | 4.15179 | 4.15171 | 4.15171 | 0.00008 | 56.2 |
| 2 | (1 1 0) | 30.422 | 30.426 | 30.426 | -0.004 | 2.93576 | 2.93541 | 2.93541 | 0.00035 | 100.0 |
| 3 | (1 1 1) | 37.489 | 37.499 | 37.499 | -0.010 | 2.39704 | 2.39643 | 2.39643 | 0.00061 | 61.1 |
| 4 | (2 0 0) | 43.562 | 43.569 | 43.569 | -0.007 | 2.07590 | 2.07558 | 2.07558 | 0.00032 | 30.1 |
| 5 | (2 1 0) | 49.021 | 49.028 | 49.028 | -0.007 | 1.85674 | 1.85648 | 1.85648 | 0.00026 | 79.7 |
| 6 | (2 1 1) | 54.059 | 54.062 | 54.062 | -0.003 | 1.69496 | 1.69488 | 1.69488 | 0.00009 | 37.8 |
| 7 | (2 2 0) | 63.304 | 63.307 | 63.307 | -0.003 | 1.46788 | 1.46782 | 1.46782 | 0.00006 | 16.5 |
| 8 | (2 2 1) | 67.640 | 67.642 | 67.642 | -0.002 | 1.38393 | 1.38390 | 1.38390 | 0.00003 | 49.2 |
| 9 | (3 1 0) | 71.846 | 71.846 | 71.846 | -0.001 | 1.31291 | 1.31290 | 1.31290 | 0.00001 | 32.8 |
| 10 | (3 1 1) | 75.952 | 75.950 | 75.950 | 0.002 | 1.25181 | 1.25184 | 1.25184 | -0.00003 | 19.7 |
| 11 | (2 2 2) | 79.986 | 79.987 | 79.987 | -0.001 | 1.19852 | 1.19851 | 1.19851 | 0.00001 | 4.1 |
| 12 | (3 2 0) | 83.970 | 83.966 | 83.966 | 0.004 | 1.15150 | 1.15154 | 1.15154 | -0.00004 | 12.8 |
| 13 | (3 2 1) | 87.925 | 87.920 | 87.920 | 0.005 | 1.10961 | 1.10967 | 1.10967 | -0.00005 | 23.2 |
| 14 | (4 0 0) | 95.824 | 95.809 | 95.809 | 0.015 | 1.03795 | 1.03807 | 1.03807 | -0.00013 | 3.5 |
| 15 | (3 2 2) | 99.806 | 99.793 | 99.793 | 0.014 | 1.00696 | 1.00706 | 1.00706 | -0.00010 | 19.6 |
| 16 | (3 3 0) | 103.837 | 103.821 | 103.821 | 0.016 | 0.97859 | 0.97869 | 0.97869 | -0.00011 | 16.9 |
| 17 | (3 3 1) | 107.939 | 107.913 | 107.913 | 0.026 | 0.95249 | 0.95264 | 0.95264 | -0.00015 x | 6.6 |
| 18 | (4 2 0) | 112.138 | 112.111 | 112.111 | 0.028 | 0.92837 | 0.92852 | 0.92852 | -0.00015 x | 10.8 |
| 19 | (4 2 1) | 116.467 | 116.431 | 116.431 | 0.037 | 0.90600 | 0.90618 | 0.90618 | -0.00018 x | 19.2 |
| 20 | (3 3 2) | 120.966 | 120.921 | 120.921 | 0.045 | 0.88517 | 0.88536 | 0.88536 | -0.00020 x | 9.1 |
| 21 | (4 2 2) | 130.709 | 130.629 | 130.629 | 0.080 | 0.84748 | 0.84775 | 0.84775 | -0.00027 x | 5.5 |
| 22 | (4 3 0) | 136.142 | 136.035 | 136.035 | 0.107 | 0.83036 | 0.83067 | 0.83067 | -0.00031 x | 7.2 |
| 23 | (5 1 0) | 142.176 | 142.026 | 142.026 | 0.150 | 0.81423 | 0.81460 | 0.81460 | -0.00037 x | 23.7 |
| 24 | (3 3 3) | 149.176 | 148.945 | 148.945 | 0.230 | 0.79901 | 0.79946 | 0.79946 | -0.00044 x | 12.5 |




Fig. B2. θ - θ unit: Diffraction pattern of LaB₆ with the PDF card # 34-427 superimposed.

Table B.III. Lattice parameter refinement of LaB₆ data in Fig. B3 taken on the 4-axis unit after alignment.

| [D07914.raw] Lab6 on a zero bckgd plate | | | | | | | | | Cell Refinement Report | |
|---|---------|---------|---------|---------|--------|---------|---------|---------|------------------------|-------|
| SCAN: 15.0/155.0/0.02/0.2(sec), Cu, I(max)=4748, 06/24/01 09:11 | | | | | | | | | | |
| PEAK: 25-pts/Parabolic Filter, Threshold=1.0, Cutoff=0.05%, BG=3/1.0, Peak-Top=Summit | | | | | | | | | | |
| NOTE: Intensity = Counts, 2T(0)=0.0(°), Wavelength to Compute d-Spacing = 1.54056A (Cu/K-alpha1) | | | | | | | | | | |
| Cell Type = Cubic, Pm3m (221) | | | | | | | | | | |
| Refined Cell = 4.15659(0.000018) [] | | | | | | | | | | |
| Vol= 71.81 A^3, Density(c)= (Chemical Formula =(Unknown), Z=0.0) | | | | | | | | | | |
| 2-Theta Error Window = 0.1(°), Zero Offset = 0.0(°), Displacement = 0.0(°) | | | | | | | | | | |
| ESD of Fit = 0.003(°), Delta 2-Theta = 0.0024(°), Delta d = 0.00003(A), F(19) = 413.3(19) | | | | | | | | | | |
| Intensity Weighting = Sqrt(%), 2-Theta Range = 15.0/155.0(°), (x) Outlier Rejection at 2.0 sigmas | | | | | | | | | | |
| 2T(cor) = 2T(obs) - Zero Offset - Displacement (cal=Calculated, obs=Observed, cor=Corrected) | | | | | | | | | | |
| # | (h k l) | 2T(cal) | 2T(cor) | 2T(obs) | Delta | d(cal) | d(cor) | d(obs) | Del-d | I% |
| 1 | (1 0 0) | 21.359 | 21.337 | 21.337 | 0.022 | 4.15659 | 4.16089 | 4.16089 | -0.00430 x | 55.1 |
| 2 | (1 1 0) | 30.386 | 30.377 | 30.377 | 0.009 | 2.93915 | 2.94002 | 2.94002 | -0.00086 x | 100.0 |
| 3 | (1 1 1) | 37.444 | 37.435 | 37.435 | 0.009 | 2.39981 | 2.40036 | 2.40036 | -0.00055 x | 47.9 |
| 4 | (2 0 0) | 43.509 | 43.508 | 43.508 | 0.001 | 2.07830 | 2.07835 | 2.07835 | -0.00005 | 22.7 |
| 5 | (2 1 0) | 48.960 | 48.963 | 48.963 | -0.003 | 1.85888 | 1.85879 | 1.85879 | 0.00009 | 65.7 |
| 6 | (2 1 1) | 53.992 | 53.992 | 53.992 | 0.000 | 1.69692 | 1.69692 | 1.69692 | 0.00001 | 33.2 |
| 7 | (2 2 0) | 63.222 | 63.228 | 63.228 | -0.006 | 1.46958 | 1.46945 | 1.46945 | 0.00012 | 12.1 |
| 8 | (2 2 1) | 67.552 | 67.550 | 67.550 | 0.001 | 1.38553 | 1.38556 | 1.38556 | -0.00003 | 32.4 |
| 9 | (3 1 0) | 71.750 | 71.752 | 71.752 | -0.002 | 1.31443 | 1.31439 | 1.31439 | 0.00004 | 23.3 |
| 10 | (3 1 1) | 75.849 | 75.852 | 75.852 | -0.003 | 1.25326 | 1.25321 | 1.25321 | 0.00005 | 16.8 |
| 11 | (2 2 2) | 79.875 | 79.855 | 79.855 | 0.020 | 1.19990 | 1.20015 | 1.20015 | -0.00025 x | 3.5 |
| 12 | (3 2 0) | 83.851 | 83.851 | 83.851 | 0.000 | 1.15283 | 1.15283 | 1.15283 | 0.00000 | 9.8 |
| 13 | (3 2 1) | 87.798 | 87.794 | 87.794 | 0.004 | 1.11090 | 1.11093 | 1.11093 | -0.00004 | 21.2 |
| 14 | (4 0 0) | 95.678 | 95.677 | 95.677 | 0.001 | 1.03915 | 1.03916 | 1.03916 | -0.00001 | 4.1 |
| 15 | (3 2 2) | 99.649 | 99.652 | 99.652 | -0.002 | 1.00812 | 1.00810 | 1.00810 | 0.00002 | 19.8 |
| 16 | (3 3 0) | 103.668 | 103.670 | 103.670 | -0.001 | 0.97972 | 0.97971 | 0.97971 | 0.00001 | 16.2 |
| 17 | (3 3 1) | 107.757 | 107.763 | 107.763 | -0.006 | 0.95359 | 0.95355 | 0.95355 | 0.00004 | 6.6 |
| 18 | (4 2 0) | 111.942 | 111.940 | 111.940 | 0.002 | 0.92944 | 0.92945 | 0.92945 | -0.00001 | 10.5 |
| 19 | (4 2 1) | 116.254 | 116.256 | 116.256 | -0.002 | 0.90704 | 0.90703 | 0.90703 | 0.00001 | 21.6 |
| 20 | (3 3 2) | 120.733 | 120.733 | 120.733 | 0.000 | 0.88619 | 0.88619 | 0.88619 | 0.00000 | 9.5 |
| 21 | (4 2 2) | 130.421 | 130.418 | 130.418 | 0.004 | 0.84846 | 0.84847 | 0.84847 | -0.00001 | 7.7 |
| 22 | (4 3 0) | 135.815 | 135.822 | 135.822 | -0.008 | 0.83132 | 0.83130 | 0.83130 | 0.00002 x | 9.5 |
| 23 | (5 1 0) | 141.792 | 141.790 | 141.790 | 0.002 | 0.81517 | 0.81518 | 0.81518 | -0.00001 | 36.1 |
| 24 | (3 3 3) | 148.700 | 148.696 | 148.696 | 0.004 | 0.79994 | 0.79994 | 0.79994 | -0.00001 | 21.9 |




Fig. B3. 4-axis unit: After alignment, diffraction pattern of LaB₆ shows peak positions are near PDF card #34-427.

Table B.IV. Lattice parameter refinement of LaB₆ data in Fig. B3 taken on the θ - θ unit after alignment.

| [X01604.raw] LaB6 std., Sollers in place | | | | | | | | Cell Refinement Report | | | |
|---|---------|---------|---------|---------|--------|---------|---------|------------------------|----------|-------|--|
| SCAN: 10.0/159.99/0.01/0.2(sec), Cu, I(max)=2072, 07/06/01 21:43 | | | | | | | | | | | |
| PEAK: 21-pts/Parabolic Filter, Threshold=1.0, Cutoff=0.05%, BG=3/1.0, Peak-Top=Summit | | | | | | | | | | | |
| NOTE: Intensity = Counts, 2T(0)=0.0(°), Wavelength to Compute d-Spacing = 1.54056Å (Cu/K-alpha1) | | | | | | | | | | | |
| Cell Type = Cubic, Pm3m (221) | | | | | | | | | | | |
| Refined Cell = 4.15686(0.00013) [Å] | | | | | | | | | | | |
| Vol= 71.83 Å ³ , Density(c)= 4.7103 (Chemical Formula =LaB6, Z=1.0) | | | | | | | | | | | |
| 2-Theta Error Window = 0.1(°), Zero Offset = 0.0(°), Displacement = 0.0(°) | | | | | | | | | | | |
| ESD of Fit = 0.0143(°), IDelta 2-Theta = 0.0114(°), IDelta dI = 0.00023(Å), F(21) = 88.0(21) | | | | | | | | | | | |
| Intensity Weighting = Sqrt(I%), 2-Theta Range = 10.0/159.99(°), (x) Outlier Rejection at 2.0 sigmas | | | | | | | | | | | |
| 2T(cor) = 2T(obs) - Zero Offset - Displacement (cal=Calculated, obs=Observed, cor=Corrected) | | | | | | | | | | | |
| # | (h k l) | 2T(cal) | 2T(cor) | 2T(obs) | Delta | d(cal) | d(cor) | d(obs) | Del-d | I% | |
| 1 | (1 0 0) | 21.358 | 21.360 | 21.360 | -0.003 | 4.15686 | 4.15637 | 4.15637 | 0.00049 | 43.1 | |
| 2 | (1 1 0) | 30.384 | 30.393 | 30.393 | -0.009 | 2.93934 | 2.93849 | 2.93849 | 0.00085 | 100.0 | |
| 3 | (1 1 1) | 37.441 | 37.453 | 37.453 | -0.012 | 2.39996 | 2.39924 | 2.39924 | 0.00072 | 43.9 | |
| 4 | (2 0 0) | 43.506 | 43.523 | 43.523 | -0.017 | 2.07843 | 2.07765 | 2.07765 | 0.00079 | 37.0 | |
| 5 | (2 1 0) | 48.957 | 48.969 | 48.969 | -0.012 | 1.85900 | 1.85857 | 1.85857 | 0.00043 | 65.0 | |
| 6 | (2 1 1) | 53.988 | 53.999 | 53.999 | -0.010 | 1.69703 | 1.69673 | 1.69673 | 0.00030 | 39.6 | |
| 7 | (2 2 0) | 63.218 | 63.216 | 63.216 | 0.002 | 1.46967 | 1.46971 | 1.46971 | -0.00004 | 12.3 | |
| 8 | (2 2 1) | 67.547 | 67.549 | 67.549 | -0.002 | 1.38562 | 1.38558 | 1.38558 | 0.00004 | 34.7 | |
| 9 | (3 1 0) | 71.745 | 71.746 | 71.746 | -0.001 | 1.31451 | 1.31449 | 1.31449 | 0.00002 | 26.9 | |
| 10 | (3 1 1) | 75.843 | 75.844 | 75.844 | -0.001 | 1.25334 | 1.25333 | 1.25333 | 0.00001 | 14.6 | |
| 11 | (2 2 2) | 79.869 | 79.864 | 79.864 | 0.005 | 1.19998 | 1.20004 | 1.20004 | -0.00006 | 2.5 | |
| 12 | (3 2 0) | 83.844 | 83.842 | 83.842 | 0.003 | 1.15291 | 1.15294 | 1.15294 | -0.00003 | 9.8 | |
| 13 | (3 2 1) | 87.790 | 87.784 | 87.784 | 0.007 | 1.11097 | 1.11104 | 1.11104 | -0.00007 | 20.4 | |
| 14 | (4 0 0) | 95.670 | 95.653 | 95.653 | 0.016 | 1.03922 | 1.03935 | 1.03935 | -0.00013 | 4.0 | |
| 15 | (3 2 2) | 99.641 | 99.623 | 99.623 | 0.018 | 1.00819 | 1.00832 | 1.00832 | -0.00013 | 17.3 | |
| 16 | (3 3 0) | 103.659 | 103.641 | 103.641 | 0.018 | 0.97978 | 0.97990 | 0.97990 | -0.00012 | 12.6 | |
| 17 | (3 3 1) | 107.747 | 107.724 | 107.724 | 0.024 | 0.95365 | 0.95379 | 0.95379 | -0.00014 | 6.7 | |
| 18 | (4 2 0) | 111.931 | 111.908 | 111.908 | 0.023 | 0.92950 | 0.92963 | 0.92963 | -0.00013 | 8.2 | |
| 19 | (4 2 1) | 116.242 | 116.219 | 116.219 | 0.023 | 0.90710 | 0.90722 | 0.90722 | -0.00011 | 17.5 | |
| 20 | (3 3 2) | 120.720 | 120.694 | 120.694 | 0.025 | 0.88625 | 0.88636 | 0.88636 | -0.00011 | 7.5 | |
| 21 | (4 2 2) | 130.405 | 130.371 | 130.371 | 0.035 | 0.84852 | 0.84863 | 0.84863 | -0.00012 | 4.8 | |
| 22 | (4 3 0) | 135.796 | 135.788 | 135.788 | 0.008 | 0.83137 | 0.83140 | 0.83140 | -0.00002 | 4.9 | |
| 23 | (5 1 0) | 141.771 | 141.725 | 141.725 | 0.046 | 0.81523 | 0.81534 | 0.81534 | -0.00011 | 20.2 | |
| 24 | (3 3 3) | 148.673 | 148.622 | 148.622 | 0.051 | 0.79999 | 0.80009 | 0.80009 | -0.00010 | 11.4 | |




Fig. B4. θ - θ unit: After alignment, diffraction pattern of LaB₆ shows peak positions are near PDF card #34-427.

should be <20 mm at the lowest 2θ of the scan in order to keep the beam on the sample. (If several slit combinations are used, two or three scans with different slits may be required.)

3. Collect a θ - 2θ powder diffraction pattern in continuous scan mode from $\sim 15^\circ 2\theta$ (Cu $K\alpha$ radiation assumed) to the maximum possible diffraction angle of the unit at a step size of 0.02° and a total scan time of at least 4 hours.

4. Use profile-fitting software to obtain the position of each of the diffraction peaks.

5. Use the cell refinement option of your software to determine the calibration curve and save the calibration parameters.

6. Compare the calibration results with those obtained previously. If changes greater than 0.02° occur in any of the peaks, then a reason for the change must be determined before assuming that the calibration, sample mount, or instrument alignment is acceptable.

TEST RECORD

1. Record that the test was performed and indicate the average deviation from the calibration curve in the instrument maintenance and

instrument/calibration logbook. Date and sign each entry.

Appendix C. Intensity, FWHM/Resolution and X-ray Wavelength Contamination Tests

PURPOSE

The diffracted intensity and profile breadth tests are used to monitor both the performance of the X-ray tube and the alignment of the goniometer.

BACKGROUND

The intensity and spectrum of the X-ray beam produced from an X-ray tube deteriorates with normal use over time. The high-energy electrons impinging on the target cause erosion of the anode surface, producing a crater. As this depression grows, the intensity of X-rays produced is reduced by the partial absorption of the X-ray by the crater's shoulder. This results in a loss of diffracted intensity and potentially an increase in the breadth of the peak profile as the target focal area becomes larger.

In addition, the high temperature of the tungsten cathode leads to vaporization of tungsten;




Fig. C1. 4-axis unit: Diffraction pattern of (101) quartz shows reasonable peak intensity and no W or $k\beta$ lines.




Fig. C2. 4-axis and θ - θ units: The (A) mitten and (B) five fingers of quartz, respectively.

tungsten then deposits on colder surfaces inside the tube, including the target and the beryllium windows. When this occurs, tungsten L-alpha radiation is produced and emitted in the beam along with the target K-alpha lines and the Bremsstrahlung. This deposited tungsten also reduces the intensity of the target K-alpha lines by absorption of the electron beam as well as the excited X-rays and cannot be discriminated against by the electronics because of the close wavelength of the two.

PROCEDURE

1. Set generator (kV and mA) to normal power level.
2. Install typical slits and record their values.
3. For a $\theta-2\theta$ goniometer with single sample holder attachment mount the quartz plate in the standard manner. For a $\theta-\theta$ goniometer with a high temperature furnace attachment mount the heater reference strip with quartz sample attached and adjust the height of the chamber until the (101) peak is in correct two-theta position.
4. Collect data from a 5° range centered on the (101) and the (212) theoretical peak locations with a maximum step size of 0.01° and a continuous scan rate of $<1.0^\circ/\text{min}$.
5. Plot the data for the (101) peak so the peak intensity is full scale. Examine the region for evidence of contaminant radiation peaks. For copper radiation, $\text{Cu K}\beta$ occurs at 24.04° and the $\text{W L}\beta$ occurs at $25.52^\circ 2\theta$.
6. Profile fit the raw data for each of the regions and compare the intensity of the (101) peak with those recorded earlier. If the intensity drops dramatically from the previous value, or has declined to less than 50% of the tube's original values, then replacement of the X-ray tube should be considered.
7. Calculate the figure of merit of the resolution of the "five fingers of quartz" [2]. The figure of merit (FOM) is the intensity of the (212) minus the background intensity divided by the average intensity of the valleys surrounding the (212) line [i.e., trough between $(212)_{\alpha_1-\alpha_2'}$, $(212)_{\alpha_2-}$, $(203)_{\alpha_1'}$, $(301)_{\alpha_1-\alpha_2}$] minus the background intensity. *Generally, an acceptable value for the FOM should be greater than 2.*

Appendix D. Detector Linearity Calibration (Dead Time Correction)

PURPOSE

The purpose of Detector Linearity Calibration is to measure the detector linearity at various count rates, and determine calibration param-

eters to correct for the intensity data that exceeds the detector linearity limit. The lack of dead-time correction can lead to improper intensity and FWHM values for intense peaks, which partially saturate the detector.

BACKGROUND

At high count rates, some X-ray detectors (e.g., liq. N_2 cooled HPGe and PSD detectors) and associated electronics cannot process each incident photon fast enough to prevent overlap of signals. Missing these counts leads to errors in intensity readings. In other words, at very high count rates, the number of events recorded by the detector is lower than the true number of incident events. Measuring the dead-time parameter permits the system to correct the observed count rate for these effects, up to a maximum observed count rate established during the calibration.

CAUTION!

The measurement of the dead-time should only be performed by the custodian of the X-ray system. The power settings (kV times mA) for the particular X-ray tube should not exceed the recommended maximum settings provided by the tube manufacturer.

PROCEDURE

- The dead-time correction should be determined for each detector on an annual basis, or when experimental conditions have changed (i.e. repair or replacement of detector or detector electronics, etc.).
- The dead-time correction should be determined as per manufacturer software and on-line help if available. If unavailable, the dead-time correction function should be determined manually and recorded.

For Manufacturer Software Equipped Systems:

- The maximum allowable potential (kV) setting for the particular X-ray tube is set, and the current (mA) setting is incremented per software window instructions.

For Systems without Software Packages:

- The dead time, t_d , of the detector and associated electronics is determined by measuring count rate (i.e., measured intensity, I_m) as a function of generator current for the (101) and (202) reflections of quartz.
- At least four scans of these reflections should be made at four "low" generator currents, which result in low count rates (typically $<10^4$ cps).
- These data should be in the linear portion of




Fig. D1. Measured intensity, I_m , as a function of generator current for the (101) reflection of quartz (linear portion of the curve).




Fig. D3. Measured intensity as a function of generator current for the (202) reflection of quartz.




Fig. D2. Measured and true intensities as a function of generator current for the (101) reflection of quartz.




Fig. D4. $\ln(I_m/I_t)$ as a function of I_t for the (101) reflection of quartz (non-linear portion of the curve).

the count rate versus generator current curve.

- These data are then fitted with straight lines in order to predict the true count rates, I_t , at larger generator currents, where dead time errors can occur for high intensity reflections (see Fig. D1).
- Both reflections should also be measured at four “high” generator currents, such that the count rates of the (101) and (202) reflections were non-linear and linear, respectively, with generator current (see Figs. D2 and D3, respectively). This “linear” count rate for the (202) reflection will demonstrate that the detector can be linear at higher generator currents.
- The dead time can be estimated using the following relation:

$$I_m/I_t = \exp(-I_t t_d),$$

where I_t is the true intensity.

- The dead-time is determined from the slope of Fig. D4 and in our example was found to be $6.9 \mu\text{s}$. The data is then corrected using [2]:

$$I_{\text{corr}} = I_m / (1 - I_m t_d)$$

Example Measurement

- The example intensities were measured using the dead-time correction program in DMSNT as well as independently. The results from each data set were similar. In Table D.I, the true intensity originates from the linear extrapolation of the data as shown in Fig. D2. The percent fit can be improved if the True intensity is plotted as a function of Measured and fit with either a power law or polynomial function.

Table D.I. The raw and predicted dead-time data.

| Generator Current (mA) | True Intensity (cps) | Measured Intensity (cps) | Percent Difference¥ | "Manual" Correction* | Percent Fit¥¥ |
|---------------------------|----------------------------|--------------------------------|------------------------|-------------------------|------------------|
| 400 | 63341 | 43793 | 31 | 60115 | 5 |
| 360 | 57035 | 41601 | 27 | 56060 | 2 |
| 320 | 50729 | 38625 | 24 | 50787 | 0 |
| 280 | 44422 | 35556 | 20 | 45611 | -3 |
| 240 | 38116 | 31911 | 16 | 39782 | -4 |
| 200 | 31809 | 27768 | 13 | 33543 | -5 |
| 160 | 25503 | 23320 | 9 | 27262 | -7 |
| 120 | 19197 | 18210 | 5 | 20528 | -7 |
| 80 | 12890 | 12713 | 1 | 13801 | -7 |
| 60 | 9737 | 9701 | 0 | 10322 | -6 |
| 50 | 8160 | 8187 | 0 | 8625 | -6 |
| 40 | 6584 | 6602 | 0 | 6884 | -5 |
| 30 | 5007 | 5035 | -1 | 5197 | -4 |
| 20 | 3431 | 3394 | 1 | 3467 | -1 |

¥ Percent Difference = $(I_t - I_m) / I_t$

* Manual correction = $I_m / (1 - I_m t_d)$

¥¥ Percent Fit = $(I_t - I_{\text{predicted}}) / I_t$

Table D.I. (cont'd).

| Power Law Fit** | Percent Fit¥¥ | Polynomial Fit*** | Percent Fit¥¥ |
|--------------------|------------------|----------------------|------------------|
| 58412 | 8 | 62294 | 2 |
| 54840 | 4 | 57441 | -1 |
| 50060 | 1 | 51183 | -1 |
| 45218 | -2 | 45131 | -2 |
| 39590 | -4 | 38469 | -1 |
| 33371 | -5 | 31594 | 1 |
| 26928 | -6 | 25036 | 2 |
| 19871 | -4 | 18555 | 3 |
| 12777 | 1 | 12840 | 0 |
| 9165 | 6 | 10262 | -5 |

** Power Law Fit = $A I_m^B = 0.11565 I_m^{1.2288}$

*** Polynomial Fit = $A + B I_m + C I_m^2 = 1072 + 0.6515 I_m + 1.75 \times 10^{-5} I_m^2$