

X-ray Seamless Pixel Array Detector

XSPA-200 ER

—High energy resolution detector
for a benchtop X-ray diffractometer—



Abstract

In X-ray diffraction measurements using a Cu source, transition metals in the sample—for example, batteries and steel materials—generate fluorescent X-rays. These fluorescent X-rays raise background intensities in the measured data, making it difficult to detect peaks derived from trace crystalline phases. The new “XSPA-200 ER” detector, which can be mounted on a benchtop X-ray diffractometer, has high energy resolution, enabling measurements with low background intensities.

1. Introduction

The XSPA-200 ER X-ray Seamless Pixel Array detector is a multidimensional pixel detector for the MiniFlex benchtop X-ray diffractometer. The most significant difference between the XSPA-200 ER and the conventional detector for MiniFlex is its energy resolution. The XSPA-200 ER is a compact detector with the same high energy resolution as the “XSPA-400 ER”⁽¹⁾, which can be installed in the SmartLab automated multipurpose X-ray diffractometer. X-ray diffraction (XRD) patterns measured by XSPA-200 ER have low background intensities because the detector discriminates fluorescent X-rays generated from the sample. The XSPA-200 ER enables the MiniFlex to obtain higher quality measurement data. This paper introduces the features of the XSPA-200 ER with actual examples.

2. Features of XSPA-200 ER

2.1. Low background measurement

CuK α radiation is generally used for XRD. When transition metals such as Cr, Mn, Fe, and Co are contained in a sample, fluorescent X-rays derived from these transition metals are generated by CuK α radiation. The energy of these X-rays is known to be close to that of the CuK α X-rays. Figure 1 shows the energy distribution of three different detectors calculated as

a normal distribution using the energy resolution for CuK α X-rays, and the fluorescent X-ray energy of each element. A narrower energy distribution indicates higher energy resolution. If the energy resolution of the detector is low, the background intensity of the diffraction pattern increases because the CuK α line cannot be discriminated from the fluorescent X-rays of transition metals. The XSPA-200 ER has a narrow energy distribution, and its energy resolution is high enough to discriminate between the CuK α line and transition metal fluorescent X-rays. Therefore, the XRD patterns obtained by the XSPA-200 ER have low background intensity.

Figure 2 shows the XRD patterns of lithium nickel cobalt manganese oxide (Li(Ni_{1/3}Co_{1/3}Mn_{1/3})O₂: NCM111), a positive electrode material for lithium ion batteries, as an example of measurement. The XRD patterns were obtained by the XSPA-200 ER and a conventional detector, D/teX Ultra2. The XRD pattern by the XSPA-200 ER showed low background intensities compared to the results from the D/teX Ultra2. The PB (peak to background) ratio calculated from the strongest peak of NCM111 was 19 with the D/teX Ultra2 and 54 with the XSPA-200 ER, resulting in an improvement of about 2.8 times. As a result, the peak of trace lithium carbonate was clearly observed by the XSPA-200 ER.




Fig. 1. The energy distribution of three detectors calculated using the energy resolution for the $\text{CuK}\alpha$ line, and the fluorescent X-rays energy of each element.




Fig. 2. Comparison of XRD patterns of NCM111 measured by two different detectors.

2.2. High reduction of $\text{K}\beta$ line

Characteristic X-rays emitted from an X-ray tube include both $\text{K}\alpha$ and $\text{K}\beta$ lines. Generally, XRD measurements use only the $\text{K}\alpha$ line for analysis. Therefore, the $\text{K}\beta$ line should be removed, leaving only the $\text{K}\alpha$ line. For example, a metal filter or a crystal monochromator is commonly used in the Bragg-Brentano parafocusing geometry⁽²⁾. Although these methods are effective, it is inevitable that the $\text{K}\alpha$ line intensity decreases when removing the $\text{K}\beta$ line. XSPA-200 ER with its high energy resolution can reduce the intensity of $\text{K}\beta$ without attenuating the $\text{K}\alpha$ line. By selecting filterless measurement, the XSPA-200 ER enables measurement with high intensity.




Fig. 3. 2D measurement data of steel material.

2.3. Multidimensional measurements of XSPA-200 ER

The XSPA-200 ER can be used as a 0D, 1D, and 2D detector. Alignment of the instrument using 0D measurement, data collection of high-quality XRD patterns by 1D measurement, and observation of Debye rings by 2D measurement can all be performed without replacing the detector. Figure 3 shows the results of a 2D measurement of a steel sample using the XSPA-200 ER. Low background data were obtained even for this steel material, which generates fluorescent X-rays derived from iron. Because the Debye rings are spotty, the sample contains coarse particles. The 2D measurement by the XSPA-200 ER enables low-background measurement and provides information on the sample condition, such as the presence of coarse particles and preferred orientation.

Table 1. XSPA-200 ER specifications.

Counting method	Direct detecting photon counting
Sensor	Silicon
Pixel size	75 μm \times 75 μm
Number of pixels	32,768 pixels
Detect area	9.6 mm \times 19.2 mm = 184.32 mm ²
Count rate	$> 1 \times 10^5$ cps/pixel
Correspond wavelength	CoK α , CuK α
Detect efficiency (at CuK α)	99%
Energy resolution (at CuK α)	340 eV (in fluorescent X-ray reduction mode)

2.4. Specifications

Table 1 shows the specifications of the XSPA-200 ER.

3. Conclusion

The MiniFlex benchtop X-ray diffractometer equipped with XSPA-200 ER enables you to obtain higher-quality data. As mentioned above, the most important feature of the XSPA-200 ER is its high energy resolution. The XSPA-200 ER is useful especially for analyzing samples containing transition metals such as Cr, Mn, Fe, and Co; for example, steel and battery materials.

References

- (1) *Rigaku Journal*, 39 (2023), No. 1, 23–26.
- (2) M. Omori: *Rigaku Journal*, 37 (2021), No. 1, 12–19.